



MADE IN THE USA 

People with back pain can
STILL LIVE LIFE.


Victory Medical
USA



The V-Force Razar is a cost effective treatment for low back pain that does not require the severity of V-force Flex.

Rigid anterior and posterior panels stabilize the spine, restrict range of motion and control pelvic tilt and rotation. The anterior panel and the strapping method provide hydrostatic lift to help relieve the load and stress on the spine. Elastic inserts give a conforming fit around the hips, and the fully adjustable, patented closure system does not rely on elastic, for an extended product life. The ultra-thin profile allows for comfortable, inconspicuous wear.

HOW THE RAZAR WORKS

When the strapping system is closed over the rigid anterior panel, it provides intra-abdominal muscle compression and provides a lift to hold the weight of the abdominal area. By lifting and holding the weight, pressure is relieved from lumbar discs reducing the associated pain. The posterior and anterior panels work together, giving opposing support to stabilize the lumbar region. The result is a supported and pressure-free lumbar spine held in a healthy, correct posture.

INDICATIONS

Acute and Chronic Low Back Pain, Muscle Strains, Spondylolisthesis, Activity-specific back pain, Mild post operative stabilization

FEATURES & BENEFITS

- Unloads pressure from lumbar discs to alleviate lower back pain and discomfort
- Anterior panel creates true hydrostatic lift of lower abdominal muscles
- Rigid posterior and anterior panels stabilize the lower spine and control rotation and pelvic tilt
- 11" high back panel covers L1-L5 vertebrae
- Elastic insets provide a conforming fit around the hips
- Soft, padded V-Foam lining for maximum comfort and compliance
- Ultra-thin low profile for comfortable, inconspicuous wear
- Fully adjustable, closure system does not rely on elastic for extended product life
- Provides pain relief during activity and allows for immediate participation in rehabilitation exercises

SIZING CHART

Measure around waist at navel

SM	24"-28"
MEDIUM	29"-33"
LARGE	34"-39"
X-LARGE	40"-44"
XX-LARGE	45"-49"
XXX-LARGE	50"-55"

For additional information regarding maintenance, care and/or warranty, please visit **www.victorymedicalusa.com**.

MODEL# VFR2012

APPROVED PDAC CODE L0627

2/13