



V
TRANSFORMER

People with back pain can
STILL LIVE LIFE.


Victory Medical
USA



MODEL# VT14X2013 (Short Version) MODEL# VT14XS2013 Approved PDAC Code L0639
 MODEL# VT2012 (Short Version) MODEL# VTS2012 Approved PDAC Code L0637

The V-Transformer has lumbar support that combines rigid anterior-posterior and media-lateral control with abdominal compression. The strapping system allows for independent adjustment of the upper and lower sections, making this a very easy brace to wear.

Anatomically designed around the hips and rib line for a comfortable, custom-like fit in an off-the-shelf design, the posterior extends from sacrococcygeal junction to T-9 vertebra, and the anterior extends from symphysis pubis to xyphoid...therefore capturing even more soft tissue for ultimate support.

HOW THE V-TRANSFORMER WORKS

This brace provides intra-abdominal muscle compression and provides lift to hold the weight of the abdominal area. By lifting and holding the weight, pressure is relieved from lumbar discs, reducing the associated pain.

INDICATIONS

Bulging Discs, Spinal Stenosis, Facet Syndrome, Degenerative Disc Disease, Herniated Discs, Mechanical Back Pain, Spondylolisthesis, Spondylosis, Pain in Mid Back, Post Operative Stabilization, Compression Fracture

The posterior and anterior shells work together to give opposing support to stabilize the lumbar region. The result is a supported and pressure-free lumbar spine held in a healthy, correct posture.

FEATURES & BENEFITS

- Utilizes a compressed system which maximizes rigid support
- System allows independent adjustment of upper and lower segments of brace
- Adjustments are easy and quick
- Fabric fastener makes it easy for patients to don and doff
- Extreme patient comfort due to lightweight, breathable and washable fabric

***FEATURING V-FOAM TECH**



Our V-Foam Tech is proprietary fabric and padding, designed to breathe, keeping areas cool and comfortable for better patient compliance.

SIZING CHART FOR VFORCE-FLEX	
<i>Measure around waist at navel</i>	
XXSM	22"-26"
XSM	25"-28"
SM	25"-31"
MEDIUM	31"-35"
ML.....	31"-35"
LARGE.....	35"-42"
X-LARGE	42"-46"
2XLARGE.....	45"-51"
3XLARGE.....	50"-56"
4XLARGE.....	55"-61"
5XLARGE.....	61"-66"

SIZING CHART FOR V-TRANSFORMER	
ABDOMINAL PANEL	
<i>Measure patient waist at navel & divide by 2</i>	
<i>Interior Panel Measurement: 17"</i>	
SM	28"-32"
MEDIUM	32"-42"
LARGE.....	43"-49"
X-LARGE	50"-66"

Patient must be fitted for optimal results.